



# Therapeutic Excellence **FOR FASTER RECOVERY**

LymphaTouch® — state of the art pulsating negative pressure and mechanical vibration to accelerate lymphatic drainage and to relieve swelling.



**SAFE AND  
COMFORTABLE**



**LIGHT WEIGHT  
AND EASY TO USE**



**ADJUSTABLE  
FUNCTIONS**



**MEASURABLE  
RESULTS**

**LYMPHATOUCH.COM**

# Customized treatment **WITH NEGATIVE PRESSURE**




## OPERATING PRINCIPLE

**LymphaTouch®** expands and stretches the skin and the tissue pulling anchor filaments to dilate the endothelial openings of the lymphatic capillaries.<sup>1</sup>

Simultaneously, the expansion of fascial and connective tissue structures create space for blood circulation and lymph flow.



**IMPROVES  
SECONDARY  
LYMPHEDEMA<sup>2</sup>**



**IMPROVES  
LYMPHATIC  
CIRCULATION<sup>3</sup>**



**RELIEVES FASCIA  
STRUCTURES AND  
MUSCLE SPASMS<sup>4</sup>**

## References

1. T. Iivari, R. K. Korhonen, and J. S. Jurvelin (2013) Experimental and computational analysis of soft tissue mechanical response under negative pressure in forearm, *Skin Research & Technology*, John Wiley & Sons, DOI: 10.1111/j.1600-0846.2012.00652.x, Vol 19. Issue 1, pp. e356-e365
2. Ville-Pekka Vuorinen, Jarkko Iivari, Jukka Jurvelin, Olavi Airaksinen (2013) Lymphatic therapy using negative pressure, A clinical study with the LymphaTouch device, Research Report #5320003/221, Finnish Funding Agency for Technology and Innovation
3. Vuorinen V-P, Airaksinen Olavi, (2009) A new vacuum suction device for management of lymphedema, Program & Abstract book pp. 118., 22nd International Congress of Lymphology, 21-25 September 2009, Sydney, Australia
4. Internal data on file

LymphaTouch Inc. ISO 13485:2003 certified company. LymphaTouch and PhysioTouch units are CE certified, and FDA Class 1 Exempt devices. RX Only. LymphaTouch® is a registered trademark of LymphaTouch Inc. All rights reserved. v1 30/1767